



OMNIS titration vessel lid for 601406XXX titration vessels for Polytron

6.01405.040

Titration vessel lid for volumetric Karl Fischer titration with an OMNIS titrator with Polytron.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.01405.040

Qt.	Order no.	Description
1 PCS	6.1454.130	Form seal made of FPM for 61414xxx
		for KF and VA instruments



Optional accessories

Order no.	Description
-----------	-------------

6.02709.030 Screw nipple

Plastic screw nipple, suitable for 6014050X0 titration vessel lids.



6.02709.020 Septum stopper

Plastic septum stopper, suitable for 6014050X0 titration vessel lids.



6.02709.040 Stopper

Stopper for titration vessel lids 6014050X0.

