



# OMNIS 2 mL cylinder unit without accessories

6.01503.120

Intelligent 2 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Without dosing tubing and antidiffusion tip.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Scope of delivery 6.01503.120

| Qt.   | Order no.  | Description                                                                                 |
|-------|------------|---------------------------------------------------------------------------------------------|
| 1 PCS | 6.2803.010 | Grease (2 g)                                                                                |
|       |            | Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX |

