



OMNIS 50 mL cylinder unit without accessories

6.01503.250

Intelligent 50 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Without dosing tubing and antidiffusion tip.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.01503.250

Qt.	Order no.	Description
1 PCS	6.2803.010	Grease (2 g)
		Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX

