



Cylinder unit OMNIS special 10 mL, without accessories

6.01508.210

Intelligent cylinder unit 10 mL for one OMNIS Titrator, Titration Module or Dosing Module. This cylinder unit is especially recommended for the following solutions:

- Aqueous alkali solutions
- Titrant 5
- Silver nitrate solutions
- Nonaqueous alkali solutions
- Permanganate solutions
- EDTA solutions

Without dosing tubing and antidiffusion tip.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.01508.210

Qt.	Order no.	Description
1 PCS	6.2803.010	Grease (2 g)
		Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



Optional accessories

Order no.	Description
-----------	-------------

6.05700.460 Consumable Kit, cylinder unit OMNIS special 10 mL

Cylinder unit special 10 mL, complete with cylinder, piston, cylinder base, valve disk, centering tube, housing, and data chip. For 6.03004.210 OMNIS cylinder unit.



6.05700.140 Consumable Kit, cylinder unit OMNIS

FEP tubing connections M6 (40 cm) and antidiffusion tip for 6.03001.xxx and 6.03004.xxx OMNIS cylinder units.

