



Antidiffusion tip

6.01543.120

Titration tip with antidiffusion valve, suitable for OMNIS Karl Fischer titrator.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.01543.120

Qt.	Order no.	Description
1 PCS	6.2726.090	Antidiffusion valve For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140, 6.1543.150, 6.1543.160, 6.1543.200

