



Bottle cap, multi-use

6.01601.000

Reusable bottle cap with 3S technology for connecting reagent bottles to an OMNIS Titrator, Titration Module or Dosing Module.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.01601.000

Qt.	Order no.	Description
1 PCS	6.02701.010	PTFE sealing ring for 601601000 (5 pcs)
		PTFE sealing ring for bottle cap, multi-use.



1 PCS

6.1819.020

Tube

For M6 threaded holes. For 6.1602.1XX bottle attachments

