



OMNIS sample rack, 9 x 200 mL, (PP)

6.02041.020

OMNIS sample rack for OMNIS Sample Robot Pick&Place, suitable for 9 sample beakers. The following sample beakers can be used: 6.1459.310.

Plastic: Polypropylene (PP)

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.02041.020

Qt.	Order no.	Description
-----	-----------	-------------

The holder for rack labels permits durable labeling of the rack with an unambiguous name.



Optional accessories

Order no.	Description	
6.1459.310	Sample beaker, disposable, plastic (PP), 200 mL, 1000 pieces Sample beaker, disposable, plastic (polypropylene, PP), 200 mL, suitable for 6.02041.020 OMNIS Software sample rack, 6.2041.370 und 6.2041.830 sample racks, 1000 pieces.	
6.01404.020	Beaker adapter for 200 mL sample beaker Beaker adapter for centering the 200 mL sample beaker in the slide of the Pick&Place modules.	
6.02064.020	Retainer for 200 mL sample beaker Retainer for fixing the 200 mL sample beaker in place in the gripper fingers of the OMNIS Sample Robot Pick&Place.	