




OMNIS sample rack, 16 x 120 mL, (PP)

6.02041.030

OMNIS sample rack for OMNIS Sample Robot Pick&Place, suitable for 16 sample beakers. The following sample beakers can be used: 6.01400.200, 6.01400.300, 6.01400.303.

Plastic: Polypropylene (PP)

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.02041.030

Qt.	Order no.	Description
-----	-----------	-------------

The holder for rack labels permits durable labeling of the rack with an unambiguous name.



Optional accessories

Order no.	Description	
6.01400.200	Sample beakers, plastic (PP), 120 mL, 20 pieces Sample beakers, plastic (polypropylene, PP), 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.	
6.01404.030	Beaker adapter for OMNIS 120 mL sample beaker Beaker adapter for centering the 120 mL sample beaker in the slide of the Pick&Place modules	
6.01400.300	Sample beakers, clear glass, 120 mL, 20 pieces Sample beakers, clear glass, 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.	
6.01400.303	Sample beakers, amber glass, 120 mL, 20 pieces Sample beakers, amber glass, 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces	