




OMNIS sample rack, 25 x 75 mL, (PP)

6.02041.040

OMNIS sample rack for OMNIS Sample Robot Pick&Place, suitable for 25 sample beakers. The following sample beakers can be used: 6.01402.000, 6.01402.003, 6.1459.400.

Plastic: Polypropylene (PP)

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.02041.040

Qt.	Order no.	Description
-----	-----------	-------------

The holder for rack labels permits durable labeling of the rack with an unambiguous name.



Optional accessories

Order no.	Description	
6.01404.040	Beaker adapter for OMNIS 75 mL sample beaker Beaker adapter for centering the 75 mL sample beaker in the slide of the Pick&Place modules	
6.01402.000	Sample beakers, clear glass, 75 mL, 30 pieces Sample beakers, clear glass, 75 mL, suitable for 602041040 OMNIS sample rack, 30 pieces	
6.01402.003	Sample beakers, amber glass, 75 mL, 30 pieces Sample beakers, amber glass, 75 mL, suitable for 602041040 OMNIS sample rack, 30 pieces	
6.1459.400	Sample beaker, plastic (PE), 75 mL, 144 pieces Sample beaker, plastic (polyethylene, PE) 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, and 6.2041.800, suitable for 6.02041.040 OMNIS Software sample rack, 144 pieces.	