



OMNIS sample rack, 9 x 150 mL, (PP)

6.02041.050

OMNIS sample rack for OMNIS Sample Robot Pick&Place, suitable for 9 sample beakers. The following sample beakers can be used: 150 mL beaker (tall shape and without spout).

Plastic: Polypropylene (PP)

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.02041.050

Qt.	Order no.	Description
-----	-----------	-------------

The holder for rack labels permits durable labeling of the rack with an unambiguous name.



Optional accessories

Order no.	Description
6.01404.050	Beaker adapter for 150 mL sample beaker
	Beaker adapter for centering the 150 mL sample beaker in the slide of the Pick&Place modules.

