



OMNIS sample rack 77 x 2R

6.02041.140

OMNIS sample rack for OMNIS Sample Robot Oven, suitable for 77 2R sample vials.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.02041.140

Qt.	Order no.	Description
1 PCS	6.02042.000	Holder for rack label (3x)

The holder for rack labels permits durable labeling of the rack with an unambiguous name.

