



Gripper fingers 48-64 mm

6.02601.030

Gripper fingers for OMNIS Pick&Place Sample Robot for gripping sample beakers with an outer diameter of 48 - 64 mm.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.05700.250	Consumable Kit OMNIS Gripper for 6.02601.030
	Four rubber gripper fingers for replacement on the gripper finger of the OMNIS Sample Robot.

