



Septum stopper

6.02709.020

Plastic septum stopper, suitable for 6014050X0 titration vessel lids.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.02709.020

Qt.	Order no.	Description
1 PCS	6.1448.010	Septum 12 mm, 5 pieces
		Set of 5 items.



1 PCS

6.2730.010

Screw nipple

Made of PP. For 6.1448.010 Septum or for mounting 6.1403.040 Drying tube or 6.0338.100 Electrode in 6.1414.030 and 6.1414.040 KF Titration vessel lids. Used for KF titrations.

