




KF Dis-Cover lid for OMNIS 75 mL sample beaker, 25 pieces

6.02710.000

Lid with seals for volumetric Karl Fischer Titration in 6.01402.000 and 6.01402.003 sample beakers in the OMNIS Sample Robot Pick&Place system.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.02710.000

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.05700.160

Consumable Kit seals for the 6.02710.000 KF Dis-
Cover lid

Replacement seals for 6.02710.000 lid, 25 pieces.

