




KF Dis-Cover lid for OMNIS 120 mL sample beaker, 16 pieces

6.02710.010

Lid with seals for volumetric Karl Fischer Titration in 6.01400.200, 6.01400.300 and 6.01400.303 sample beakers in the OMNIS Sample Robot Pick&Place system.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.02710.010


Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.05700.170

Consumable Kit seals for the 6.02710.010 KF Dis-
Cover lid

Replacement seals for 6.02710.010 lid, 16 pieces



Optional accessories

Order no.	Description	
6.01400.200	Sample beakers, plastic (PP), 120 mL, 20 pieces Sample beakers, plastic (polypropylene, PP), 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.	
6.01400.300	Sample beakers, clear glass, 120 mL, 20 pieces Sample beakers, clear glass, 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.	
6.01400.303	Sample beakers, amber glass, 120 mL, 20 pieces Sample beakers, amber glass, 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces	