



Ultrafiltration cell 2

6.02729.110

For the Inline ultrafiltration of samples together with the 858 Professional Sample Processor / 919 IC Autosampler plus. As an alternative to the time-consuming manual filtration with a syringe filter.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description
6.2714.020	Filtration membrane reg. cellulose, 0.2 µm, 50 pieces
	Nominal pore size: 0.2 µm. Set of 50 pieces. For inline filtration.



6.2714.030 Dialysis membrane (polyamide)
For inline dialysis and filtration



6.2744.000 Pressure screw PVDF
With UNF 10/32 connection. For IC instruments, inline dialysis
(connection of PTFE capillaries)



6.2831.010 Plastic tweezers
For dialysis and filtration membranes



6.02057.120 Ultrafiltration cell 2 holder, stand-alone
Holder for the 6.02729.110 ultrafiltration cell 2



6.02057.030 Ultrafiltration cell 2 holder

Holder for mounting the 6.02729.110 ultrafiltration cell 2 on the 858 Professional Sample Processor / 919 IC Autosampler plus.



6.2714.050 Filtration membrane reg. cellulose, 0.45 µm, 50 pieces

Nominal pore size: 0.45 µm. Set of 50 pieces. For inline filtration.

