

Accessory kit: Bottle OMNIS IC, PP, 2 L

6.05000.010

OMNIS IC accessory set to install the reagent bottle onto the Reagent Organizer / Reagent Organizer Liquid Control.


Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.


These lists may be subject to change.

Scope of delivery 6.05000.010


Order Number	Description	Quantity
6.01602.100	Flaschenaufsatz GL 45, 2 x M6, UNF 10/32	1 PCS



6.01608.000	Bottle, plastic (PP), 2 L, GL 45	1 PCS
-------------	----------------------------------	-------




6.01834.010	FEP-Schlauch M6 mit Filter, 0.27 m	1 PCS
-------------	------------------------------------	-------



6.02050.000

Flaschenhalter, RO


1 PCS



6.02050.010

Flaschenbefestigung RO


1 PCS



Optional accessories

6.05801.040

**Volumensensor, ROLC-
Kit**



6.02061.000

Reagenzorganisierer

