

Zubehör-Kit: OMNIS IC AS

6.05000.050

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.


These lists may be subject to change.

Scope of delivery 6.05000.050


Order Number	Description	Quantity
6.01816.010	Silicone tubing, 8/12 mm, 2 m Drainage tubing	2 PCS



6.02040.000	AS Tray	1 PCS
-------------	---------	-------




6.02041.510	OMNIS Probenrack 40 x 5 mL	2 PCS
-------------	----------------------------	-------



6.02061.120

Tropfensammler, Extension


1 PCS



6.02743.027

Vial PS, 5 mL, 200 Stk.

1 PCS



6.02743.030

Bottle, PP, 50 mL, 12 pieces

1 PCS



The 50 mL bottles are made of polypropylene (PP) plastic and can be used for samples, standards, check standards and rinsing solutions. Production and packaging both take place under sterile and clean room conditions. This guarantees minimum contamination. The sample vessels are subjected to random checks for leachable anions and cations.

6.02743.087

Kappe PP, AD 11.5 mm, 200 Stk.

1 PCS

