



WSM accessory kit (2T/4P)

6.05001.210

This accessory kit contains all of the components required for the tubing setup for 2 rinsing and 2 aspiration pumps on the pump module of the OMNIS Workstation Module (2T/4P). In addition to preassembled tubing with Luer connectors, each of these also includes 4 press clamps and pump tubing. A tool for simplifying attachment of the Luer connectors is also included in the set.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.05001.210

Qt.	Order no.	Description
-----	-----------	-------------

4 PCS

6.01801.020

Pharm-A-line pump tubing to the peristaltic pump module

The pump tubing comes pre-equipped with adapters that permit direct installation in the peristaltic pump of the pump module. The Pharm-A-line pump tubing is used for the transport of both aqueous and nonaqueous solutions and must be operated in connection with the white press clamp in the pump module.



1 PCS

6.01805.100

PTFE tubing with Luer 60 cm

Tubing for the connection between peristaltic pump 1 or 3 in the pump module with the distributor (rinsing) on the Pick&Place module.



3 PCS

6.01805.110

PTFE tubing with Luer 80 cm

Tubing for the connection between peristaltic pump 2 or 4 in the pump module with the distributor (aspiration) on the Pick&Place module.



4 PCS

6.01805.120

PTFE tubing with Luer 2.5 m

Tubing for the connection between the canister and the peristaltic pump in the pump module.



1 PCS

6.02621.030

Luer tool

This tool simplifies the attachment and disassembly of Luer lock connectors. It can be fitted on the Luer lock connector, where its shape facilitates the attachment of the tubing equipped with Luer connectors.



4 PCS

6.02703.000

Press clamp, white, for "Peristaltic" pump module

The press clamp is used to fix and secure the pump tubing in the "Peristaltic" pump module.

