



## NH<sub>3</sub>-selective gas membrane electrode (low conc.)

6.0506.100

Combined ammonia-selective gas membrane electrode for low concentrations  $10^{-2}$  mol/L, e.g.:

- in drinking water
- in boiler feed water
- in mineral water

If necessary, the complete membrane module can be easily exchanged. They are individually tested and certified.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Scope of delivery 6.0506.100

| Qt. | Order no. | Description |
|-----|-----------|-------------|
|-----|-----------|-------------|

100 Stk

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



---

100 Stk

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



## Optional accessories

| Order no.  | Description                                                                                                                                            |                                                                                       |
|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| 6.1255.000 | Replacement Kit for NH <sub>3</sub> -selective Electrode 6.0506.100<br><br>Consisting of: replacement module (3 pieces), measuring electrolyte (50 mL) |    |
| 6.2104.020 | Electrode cable / 1 m / F<br><br>For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).                              |    |
| 6.2316.030 | Inner filling solution to 6.0506.100 electrode<br><br>Measuring electrolyte for NH <sub>3</sub> -selective electrode (50 mL)                           |  |