

OMNIS IC-Ausr.: D-INJ, DD

6.05331.010


Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.


These lists may be subject to change.

Scope of delivery 6.05331.010


Order Number	Description	Quantity
6.01562.100	Transferkapillare 2 mL, PEEK	1 PCS



6.01805.400	FEP-Schlauch M6, 1/8" x 1.59 mm, 3 m	1 PCS
-------------	--------------------------------------	-------




6.02744.010	flanschlose Verschraubung M6, 1/8", 2 Stk.	1 PCS
-------------	--	-------



6.05000.010

Accessory kit: Bottle OMNIS IC, PP, 2 L

1 PCS



OMNIS IC accessory set to install the reagent bottle onto the Reagent Organizer / Reagent Organizer Liquid Control.

6.2744.070

Pressure screw short

1 PCS



Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries