


OMNIS IC-Ausr.: MIDT-UF, PP/DD


Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery

Order Number	Description	Quantity
6.02744.010	flanschlose Verschraubung M6, 1/8", 2 Stk.	4 PCS




6.01815.000	Spiralband, 3-15 mm, 0.5 m	1 PCS
-------------	----------------------------	-------



6.05000.010	Accessory kit: Bottle OMNIS IC, PP, 2 L	1 PCS
-------------	--	-------

OMNIS IC accessory set to install the reagent bottle onto the Reagent Organizer / Reagent Organizer Liquid Control.



6.2744.070

Pressure screw short

2 PCS

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



6.2744.320

PEEK coupling M6 / UNF 10/32

2 PCS

Connecting transfer tubing (M6) with capillaries (UNF 10/32).



6.2621.030

Hexagon key 4 mm

1 PCS



6.1831.060

PEEK capillary 0.5 mm i.d., 1 m

1 PCS

For inline dialysis and for Sample Changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Filtration Sample Processor, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit



6.1831.050


PEEK capillary 0.5 mm i.d., 40 cm

1 PCS

For inline dialysis and for sample changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit, Compact VA Autosampler

6.2744.034


Coupling olive/UNF 10/32, 2x

1 PCS

Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.

6.2744.080


M6 thread / UNF 10/32 adapter

1 PCS

With M6 outer thread. For IC instruments.

6.2744.200


Adapter UNF 10/32 outer / M6 inner


1 PCS

Connecting 1/16 in. PEEK capillaries and transfer tubing (M6)

6.01532.140

LQH-Transferschlauch, FEP, AD 1/8


1 PCS



6.01532.130

Transferschlauch 5 mL, FEP, AD 1/8


1 PCS



6.01805.400

FEP-Schlauch M6, 1/8" x 1.59 mm, 3 m


1 PCS



6.01826.020

Pumpschlauch, ID 1.0 mm, weiss - weiss


1 PCS



6.01831.220

PEEK-Kap. FILTRATE OUT - INJ, 800/0.5 mm


1 PCS



6.01831.240

PEEK Kapillare 0.75 mm i.d. / 140 cm


1 PCS



6.01831.210

PEEK Kap. PP - SAMPLE IN, 450/0.5 mm

1 PCS



6.02729.110

Ultrafiltration cell 2

1 PCS



For the Inline ultrafiltration of samples together with the 858 Professional Sample Processor / 919 IC Autosampler plus. As an alternative to the time-consuming manual filtration with a syringe filter.

6.2714.020

Filtration membrane reg. cellulose, 0.2 µm, 50 pieces 1 PCS

Nominal pore size: 0.2 µm. Set of 50 pieces. For inline filtration.



6.1831.180

PEEK Capillary 0.5 mm i.d., 3 m

1 PCS



6.2744.000

Pressure screw PVDF

1 PCS

With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)



6.2831.010

Plastic tweezers

1 PCS

For dialysis and filtration membranes



6.1903.000

Stirring bar / 8 x 3 mm

1 PCS

Stirring bar with magnetic core, PTFE covering, 8 mm long with 3 mm diameter.

