


OMNIS IC-Ausr.: D-INJ, manuell

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery

Order Number	Description	Quantity
6.02816.000	Spritze 20 mL, Luer	1 PCS
		
6.1831.060	PEEK capillary 0.5 mm i.d., 1 m For inline dialysis and for Sample Changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Filtration Sample Processor, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit	1 PCS
		
6.2744.120	Luer F pressure screw adapter Connection between pressure screw and syringe.	1 PCS
		