


OMNIS IC-Ausr.: PCR, DD


Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery


Order Number	Description	Quantity
6.02744.010	flanschlose Verschraubung M6, 1/8", 2 Stk.	2 PCS



6.02837.040	Filterrohr	2 PCS
-------------	------------	-------



6.05000.010	Accessory kit: Bottle OMNIS IC, PP, 2 L	2 PCS
-------------	--	-------



OMNIS IC accessory set to install the reagent bottle onto the Reagent Organizer / Reagent Organizer Liquid Control.

6.2744.070

Pressure screw short

1 PCS


Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



6.02027.030

Clip Befestigungsschiene, 2 x AD 1/16

2 PCS



6.2744.150

T connector UNF 10/32

1 PCS

T connector for 3 capillaries.



6.1831.030

PEEK capillary 0.75 mm i.d. / 3 m

1 PCS


For all IC components.



6.01805.400

FEP-Schlauch M6, 1/8" x 1.59 mm, 3 m


1 PCS



6.02006.120

Halter, Befestigungsschiene


1 PCS



6.02061.130

Tropfensammler, Unit

1 PCS



6.02845.100

PCR-Kartusche, PTFE, ID 0.5 mm

1 PCS

