



OMNIS IC Ausr.: EPOD

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery

Order Number	Description	Quantity
6.1621.000	Container 10 L As rinsing or waste container in automated systems.	1 PCS
		
6.02744.010	flanschlose Verschraubung M6, 1/8", 2 Stk.	3 PCS
		
6.05000.010	Accessory kit: Bottle OMNIS IC, PP, 2 L OMNIS IC accessory set to install the reagent bottle onto the Reagent Organizer / Reagent Organizer Liquid Control.	1 PCS
		

6.1805.530

FEP tubing / M6 / 2 m

1 PCS

With light and kink protection



6.2744.200

Adapter UNF 10/32 outer / M6 inner

1 PCS


Connecting 1/16 in. PEEK capillaries and transfer tubing (M6)



6.01602.110

Flaschenaufsatz S 40, 2x M6, UNF 10/32

1 PCS



6.01503.250

OMNIS 50 mL cylinder unit without accessories

1 PCS

Intelligent 50 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Without dosing tubing and antidiffusion tip.



6.01805.400

FEP-Schlauch M6, 1/8" x 1.59 mm, 3 m

2 PCS

