



Consumable Kit, cylinder unit OMNIS

6.05700.140

FEP tubing connections M6 (40 cm) and antidiffusion tip for 6.03001.xxx and 6.03004.xxx OMNIS cylinder units.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.05700.140

Qt.	Order no.	Description
1 PCS	6.1543.200	Titration tip / M6 thread
		Titration tip with antidiffusion valve.



2 PCS

6.1805.100

FEP tubing / M6 / 40 cm

With light and kink protection



1 PCS

6.2726.090

Antidiffusion valve

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140,
6.1543.150, 6.1543.160, 6.1543.200



Optional accessories

Order no.	Description
6.01508.220	Cylinder unit OMNIS special 20 mL, without accessories

Intelligent cylinder unit 20 mL for one OMNIS Titrator, Titration Module or Dosing Module. This cylinder unit is especially recommended for the following solutions:

- Aqueous alkaline solutions
- Titrant 5
- Silver nitrate solutions
- Nonaqueous alkaline solutions
- Permanganate solutions
- EDTA solutions

Without dosing tubing and antidiffusion tip.

