



Upgrade Kit: IC Eco

6.05800.060

Upgrade Kit for retrofitting the Eco IC with a pulsation absorber.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.05800.060

Qt.	Order no.	Description
1 PCS	6.2620.150	Pulsation absorber (metal-free) For all ion chromatographs



1 PCS

6.2744.070

Pressure screw short

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries

