



MB Retrofit Kit, SES

6.05802.010

Adaptation kit for retrofitting the optimization of the 2 mm setup for sequentially chemically suppressed ProfIC or Compact IC instruments.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.05802.010

Qt.	Order no.	Description
1 PCS	6.1825.260	PEEK sample loop 5 µL For injection valve, with 2 short PEEK pressure screws

