



Ag/AgCl reference system

6.0724.140

Silver / silver chloride reference system with SGJ 14/15 for the modular reference system in combination with the reference electrolyte vessel made of PTFE (6.1240.000) or glass (6.1240.020).

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description
6.1240.000	Electrolyte vessel / diameter 3 mm With PTFE capillary and ceramic diaphragm (3 mm)



6.1240.020 Electrolyte vessel, diameter 1 mm
With ceramic diaphragm (D type).



6.2308.020 Electrolyte 3 mol/L KCl (250 mL)
Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)

