



LL-Ag/AgCl reference electrode for CVS

6.0728.130

Inner reference system with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$. For use with electrolyte vessel made of plastic 6.1245.010.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.1245.010	Electrolyte vessel made of plastic
	Electrolyte vessel made of plastic with ceramic diaphragm for double junction reference electrode.



6.2112.130

Adapter plug contact 5.5 mm / screw contact, gold-plated

For VA electrodes

