



Complete LL-Ag/AgCl(gel) reference electrode

6.0730.100

Maintenance-free Double-Junction Ag/AgCl reference electrode for voltammetry. Made of glass with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$ (gel), for pH 0 - 8. Suitable for use in continuous operation in automated systems and online instruments.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.0730.100

Qt.	Order no.	Description
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm, polyethylene.



Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



Optional accessories

Order no.	Description
6.2112.130	Adapter plug contact 5.5 mm / screw contact, gold-plated For VA electrodes

