



Ag/AgCl reference electrode (length 18.5 cm)

6.0736.110

Silver / silver chloride reference electrode, installation length 16.8 cm. This reference electrode is well suited for automated applications. The flexible ground-joint diaphragm, which is insensitive to contamination, can be replaced at any time. The recommended reference electrolyte is $c(\text{KCl}) = 3 \text{ mol/L}$. This sensor is supplied without electrolyte filling.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.2106.020

Strand / 1 m / 2 x B

Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B



6.2106.060

Strand / 2 m / 2 x B

Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B



6.2308.020

Electrolyte 3 mol/L KCl (250 mL)

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)



6.2310.010

Electrolyte KNO_3 -1 mol/L 250 mL

Electrolyte solution KNO_3 1 M (reference electrolyte for combined silver electrodes)



6.2320.000

TEABr 0.4 mol/L in ethylene glycol (250 mL)

Electrolyte solution TEABr (tetraethylammonium bromide in ethylene glycol), $c(\text{TEABr}) = 0.4 \text{ mol/L}$

