



High Throughput Lab-grade Raman Probe (532 nm), 5m

607510090

High-throughput lab-grade fiber optic Raman probe assembly with trigger button. 5 meter fiber length.

Compatible with i-Raman NxG 532.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 607510090


Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

607510430

Distance Regulator Cap

Distance regulator cap for standard lab-grade Raman probe shaft.



1 PCS

607510120

Standard Shaft for Lab-Grade Probe (532 nm)

Replacement shaft for 532 nm excitation lab-grade Raman probe assembly. Working distance 5.4 mm.

