



Lab-Grade Raman Probe (785 nm), 5 m

607510460

Lab-grade fiber optic Raman probe assembly with trigger button. 5 meter fiber length. Compatible with i-Raman Plus 785S and 785HR models.

Below, the accessories are grouped into Scope of delivery and Optional accessories.


Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 607510460

Qt.	Order no.	Description
1 PCS	607510430	Distance Regulator Cap

Distance regulator cap for standard lab-grade Raman probe shaft.



Replacement shaft for 785 nm excitation lab-grade Raman probe assembly. Working distance 5.5 mm.

