



Platinum electrode tip

6.09395.024

Electrode tip made of platinum (Pt) for Autolab RDE, 3 mm electrode disc diameter, M4 thread connection. Temperature range: 0 - 40°C.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.09395.024

Qt.	Order no.	Description
1 PCS	6.2802.000	Polishing set for solid-state electrodes
		Polishing set with 1 polishing cloth and about 2 g aluminum oxide powder (grain size 0.3 µm)

