



# Metrosep A Supp 5 S-Guard/4.0

6.1006.540

The Metrosep A Supp 5 S-Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contaminations from the sample or eluent.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Scope of delivery 6.1006.540

| Qt.   | Order no.  | Description                                                     |
|-------|------------|-----------------------------------------------------------------|
| 2 PCS | 6.2744.060 | Threaded stopper                                                |
|       |            | For UNF 10/32. Stopper for IC, e.g. for the sealing of columns. |

