



## Metrosep C 4 - 150/2.0

6.1050.220

The Metrosep C 4 - 150/2.0 "Microbore" cation separation column is the universal standard column in cation analysis of "Microbore" separating columns. High separating efficiency in a short time. The Metrosep C 4 - 150/2.0 is the ideal separation column for the analysis of alkaline and earth alkaline metals in aqueous media. With its low eluent flow, this column is particularly suitable for IC-MS coupling.

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 6.1050.220

| Qt. | Order no. | Description |
|-----|-----------|-------------|
|-----|-----------|-------------|

2 PCS

6.2744.060

Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



## Optional accessories

| Order no.  | Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |                                                                                      |
|------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| 6.1050.600 | Metrosep C 4 Guard/2.0                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |   |
|            | <p>The Metrosep C 4 Guard/2.0 contains the C-4 column material and is used to protect all Metrosep cation columns with 2 mm inner diameters which have a substrate based on silica gel. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The economical price is an additional plus. The Metrosep C 4 Guard/2.0 also works based on the "On Column Guard System" and is screwed directly onto the separation column with virtually no dead volume.</p> |                                                                                      |
| 6.1050.610 | Metrosep C 4 S-Guard/2.0                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |  |