



# Metrosep C 4 - 150/4.0

6.1050.420

The Metrosep C 4 - 150/4.0 is the universal standard column in cation analysis for accomplishing high separating efficiency in a short time. The Metrosep C 4 - 150/4.0 is the ideal separation column for the analysis of alkaline and earth alkaline metals in aqueous media.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

## Scope of delivery 6.1050.420

| Qt.   | Order no.  | Description                                                     |
|-------|------------|-----------------------------------------------------------------|
| 2 PCS | 6.2744.060 | Threaded stopper                                                |
|       |            | For UNF 10/32. Stopper for IC, e.g. for the sealing of columns. |



## Optional accessories

| Order no.  | Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |                                                                                       |
|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| 6.1050.500 | Metrosep C 4 Guard/4.0                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |    |
|            | <p>The Metrosep C 4 Guard/4.0 contains the C-4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The economical price is an additional plus. The Metrosep C 4 Guard/4.0 also works based on the "On Column Guard System" and is screwed directly onto the separation column with virtually no dead volume.</p> |                                                                                       |
| 6.1050.510 | Metrosep C 4 S-Guard/4.0                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |   |
|            | <p>The Metrosep C 4 S-Guard/4.0 contains the C-4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The economical price is an additional plus.</p>                                                                                                                                                             |                                                                                       |
| 6.1050.530 | Metrosep C 4 S-Guard - 50/4.0                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |  |
|            | <p>The Metrosep C 4 S-Guard - 50/4.0 contains the C 4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. In comparison with Standard C 4 guard columns, the Metrosep C 4 S-Guard - 50/4.0 has a greater capacity and therefore an even longer service life.</p>                                                 |                                                                                       |