



## Pt1000 temperature sensor (installation length 12.5 cm)

6.1110.100

Pt1000 temperature sensor (class B) made of glass.

This PT1000 temperature sensor is also available under the article number 6.1110.110 with an installation length of 17.8 cm.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Scope of delivery 6.1110.100

| Qt. | Order no. | Description |
|-----|-----------|-------------|
|-----|-----------|-------------|

100 Stk

6.1236.020

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm and O-ring.



# Optional accessories

| Order no.  | Description                                                                                                                                                    |                                                                                     |
|------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| 6.2104.140 | Electrode cable / 1 m / 2 x 2 mm<br><br>For connecting temperature sensors with Metrohm plug-in head G with 780, 781, 826, 827 pH/Ion meters and all Titrandos |  |
| 6.2104.150 | Electrode cable 2 m, 2 x 2 mm<br><br>For connecting temperature sensors with Metrohm plug-in head G with 780, 781, 826, 827 pH/Ion meters and all Titrandos    |  |