



Pt100 temperature sensor for stability measuring devices (30 cm)

6.1111.020

Pt100 temperature sensor (length 30.0 cm) with fixed cable and mini-DIN-plug. This glass-free temperature sensor is well suited for use with instruments for stability measurement (873 Biodiesel Rancimat and 893 Professional Biodiesel Rancimat).

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.2042.040

Clamp for temperature sensor

Holding clamp for Pt100 temperature sensor for stability measuring instruments

