



Optrode M2

6.1125.000

Optical sensor for photometric titrations offering 8 different wavelengths. The wavelength can be switched using the software (tiamo 2.5 or higher) or with a magnet. The glass shaft is completely solvent-resistant and easy to clean. For example, this space-saving sensor is suitable for:

- Non-aqueous titrations in accordance with USP or EP
- Determinations of carboxyl end groups
- TAN/TBN in accordance with ASTM D974
- Sulfate determination
- Fe, Al, Ca in cement
- Water hardness
- Chondroitin sulfate in accordance with USP

The sensor is not suitable for determinations of concentrations via measurement of color intensity (colorimetry).

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.1125.000

Qt.	Order no.	Description
100 Stk	6.1903.030	Stirring bar / 25 mm
		Stirring bar with magnetic core, PTFE covering, length 25 mm.
		
100 Stk	6.2008.040	Storage vessel
		Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.
		