



Lateral Gold Electrode Tip

6.1204.150

Metal electrode tip, gold (Au), with lateral electrode surface. Electrode disk diameter 3 mm. M3 thread. Temperature range: 0 ... 40 °C.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.2802.000	Polishing set for solid-state electrodes Polishing set with 1 polishing cloth and about 2 g aluminum oxide powder (grain size 0.3 µm)

