



Ultra Trace electrode tip

6.1204.180

Carbon electrode tip from Ultra Trace graphite. Electrode disk diameter 2 mm. M3 thread connection. Temperature range: 0 - 40 °C.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.2103.110	Contact pin M3, 2 mm For 6.1241.060 Electrode shaft and 6.1204.1XX Electrode tips. With M3 thread. Used with VA equipment sets and for ELCD applications.



6.2802.000

Polishing set for solid-state electrodes

Polishing set with 1 polishing cloth and about 2 g aluminum oxide powder (grain size 0.3 μm)



6.2802.020

Polishing set for Ultra Trace Electrode Tip

Polishing set to renew the surface of the Ultra Trace Electrode tip
6.1204.180



6.2802.030

Polishing tape to Polishing Set for Ultra Trace Electrode

Polishing tape to Polishing Set for Ultra Trace Electrode. 12 μm particle size. used in the polishing cylinder. Set of 50

