



# 1 mm platinum electrode tip for CVS

6.1204.190

Metal electrode tip made of platinum (Pt), diameter of electrode disk 1 mm  $\pm 0.02$  mm, polished surface, shaft made of glass, M3 thread connection. Concentricity error: 0.2 mm. Temperature range: 0 - 50 °C.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

## Optional accessories

| Order no.  | Description                                                                                                                                                             |
|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 6.2103.110 | <p>Contact pin M3, 2 mm</p> <p>For 6.1241.060 Electrode shaft and 6.1204.1XX Electrode tips. With M3 thread. Used with VA equipment sets and for ELCD applications.</p> |

