



Glassy carbon electrode tip

6.1204.300

Electrode tip made of glassy carbon (GC) for Autolab RDE, electrode disk diameter 3 mm, M4 thread connection. Temperature range: 0 - 40°C.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1204.300

Qt.	Order no.	Description
1 PCS	6.2802.000	Polishing set for solid-state electrodes
		Polishing set with 1 polishing cloth and about 2 g aluminum oxide powder (grain size 0.3 µm)



Optional accessories

Order no.	Description	
6.2103.120	Contact pin M4 / 2 mm For 6.1241.060 Electrode shaft and 6.1204.3XX Electrode tips. With M4 thread. Used with VA equipment sets.	
