



Electrolyte vessel / diameter 3 mm

6.1240.000

With PTFE capillary and ceramic diaphragm (3 mm)

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.0724.140	Ag/AgCl reference system Silver / silver chloride reference system with SGJ 14/15 for the modular reference system in combination with the reference electrolyte vessel made of PTFE (6.1240.000) or glass (6.1240.020).

