



Electrode holder

6.1241.020

Electrode holder for VA auxiliary electrodes with lock bolt. For use with glassy carbon rod (6.1247.000) or platinum pin (6.1247.010).

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.1247.000	Glassy carbon rod (GC) For electrode holder 6.1241.020 or 6.1241.120 for VA auxiliary electrode.

A photograph of a glassy carbon rod (GC), model 6.1247.000. It is a long, thin, dark-colored rod with a slightly tapered end.

6.1247.010

Platinum pin (Pt)

For electrode holder 6.1241.020 or 6.1241.120 for VA auxiliary electrode.

