



Electrode shaft

6.1241.040

Electrode shaft for separate metal pin electrodes. For sample changer applications.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1241.040

Qt.	Order no.	Description
1 PCS	6.2106.060	Strand / 2 m / 2 x B Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B



Optional accessories

Order no.	Description	
6.1248.000	Platinum Electrode Rod Platinum Electrode Rod for separate metal pin electrode	
6.1248.010	Silver Electrode Rod Silver Electrode Rod for separate metal pin electrode	
6.1248.030	Gold Electrode Rod Gold Electrode Rod for separate metal pin electrode	
6.1248.040	Glassy Carbon Electrode Rod Glassy Carbon Electrode Rod for separate metal pin electrode	

6.1248.050

Tungsten Electrode Rod

Tungsten Electrode Rod for separate metal rod electrode

