



Electrode holder

6.1241.120

Electrode holder for VA auxiliary electrodes with plug contact 5.5. mm. For use with glassy carbon rod (6.1247.000) or platinum pin (6.1247.010).

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.1247.000	Glassy carbon rod (GC) For electrode holder 6.1241.020 or 6.1241.120 for VA auxiliary electrode.

A photograph of a thin, dark, cylindrical glassy carbon rod (GC) shown at an angle against a light background.

6.1247.010

Platinum pin (Pt)

For electrode holder 6.1241.020 or 6.1241.120 for VA auxiliary electrode.



6.2112.130

Adapter plug contact 5.5 mm / screw contact, gold-plated

For VA electrodes

