




SPE measuring head for Professional VA instruments

6.1256.030

Measuring head for operation with *screen-printed electrodes* (SPE) or the scTRACE Gold.


Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.1256.030

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.2069.000

Park station for the measuring head of Professional
VA/CVS instruments

1 PCS

6.1829.030

PTFE tubing / M6 / 21 cm, tapered



1 PCS

6.2709.110

Stopper for the electrode opening on Professional VA
/CVS instruments

Set with 1 stopper for electrode opening (e.g. for pH electrode or Pt1000 temperature sensor), 6 stoppers for openings (auxiliary solutions and/or automatic sample addition).




1 PCS

6.1911.000

Flexible drive shaft for stirrer

10 pcs, for Professional VA instruments, including wrench.



Optional accessories

Order no.	Description	
6.1454.130	Form seal made of FPM for 61414xxx for KF and VA instruments	
6.2709.130	Adapter B-14/15 - SPE measuring head insert (VA) Adapter for use of the electrode shaft scTRACE Gold 6.1241.080 or the bismuth drop electrodes 6.0346.000 in the 884 Professional VA with SPE measuring head 6.1256.030 or in the 797 VA Computrace or 663 VA Stand with measuring head insert 6.1412.020.	